## QUICK START GUIDE BATTERY

Models: 09-0021-01, 09-0021-02, and 09-0021-03

This Quick Start Guide provides an introduction to using your Relaxis. **Please read your Relaxis Battery Operation Manual before use.** It contains important safety information, troubleshooting solutions and more.

### SET UP

1. Plug the power supply into an outlet and connect the power supply to the controller (Figure 1). Make sure to only use the provided power supply.
2. Allow the battery (within the controller) to charge while plugged into the outlet. On full charge, the battery will allow you to use Relaxis for up to two 35-minute therapy cycles without plugging the power supply into an outlet.
3. Place Relaxis where you can reach it easily from your bed.


Figure 1: Connecting the power supply to the controller

### FAMILIARIZE YOURSELF WITH RELAXIS

Learn your preferences so that when you need therapy, you can simply turn on your Relaxis, focus on the vibration, experience relief, and fall asleep.

- Following the steps below to use your Relaxis, explore which intensity feels most comfortable to you. Of course, each RLS attack may be different, and it is okay to adjust the vibration intensity to suit your needs each time you use Relaxis.
- Learn how you like to position Relaxis under or between your legs when you are lying in bed. Make sure to keep the pad flat. **(Do not bend or wrap the pad around your legs.)**
- Allow time to determine if Relaxis is improving your quality of sleep. It may take up to 4 weeks for you to see an improvement.

### USING RELAXIS

Only use Relaxis as needed. When you feel the oncoming symptoms of RLS:

1. Place Relaxis under the area of greatest discomfort and lie in the position that you find most comfortable for sleeping.
2. Turn on Relaxis by pressing the Power Button (Figure 2).
3. Use the two arrow buttons to adjust the vibration intensity to the level that provides pleasant relief (Figure 2). To change the vibration intensity rapidly, hold down on an arrow button.

*Note: An Intensity Indicator Light will either turn on or turn off when you have increased or decreased the vibration intensity eight times.*

4. Take a deep breath, and let Relaxis help you fall asleep.

*Note: During the first few minutes of use, it may seem as if the vibration “competes” with the uncomfortable RLS sensations. It could take a few minutes for you to feel pleasant relief from the vibration.*

5. After 30 minutes of vibration therapy at the intensity you set, Relaxis will automatically decrease intensity over an additional 5 minutes. At the end of the 35-minute cycle, Relaxis will shut itself off.
6. If you have not fallen asleep at the end of a 35-minute cycle, you may press the Power Button again to initiate a second cycle.
7. To stop therapy at any time, press the Power Button.

## BATTERY UPKEEP

The life and performance of the rechargeable battery depends upon how you store and care for your Relaxis.

To maximize the life of the battery:

- Store your Relaxis at room temperature. (Elevated temperatures during use and storage degrade the life and performance of the battery.)
- Recharge the battery at least once every 6 months. Otherwise, the battery may lose its ability to recharge.

*Note: If your battery loses its ability to recharge, you can still use your Relaxis while plugged into an outlet.*

## CARE AND CLEANING

Maximize the life of Relaxis through proper care and cleaning. For full details and instructions on how to care for and clean Relaxis, read your Relaxis Battery Operation Manual.

To Clean Relaxis:

1. Unplug Relaxis from the outlet and disconnect the power supply from the controller.
2. Remove the outer cotton cover and hand or machine wash (inside out) in cold water with like colors. Dry (inside out) on low or line dry.
3. Wipe the controller and pad surfaces with a slightly dampened (water only) cloth. Do not use cleaning solutions on Relaxis. (These could damage the system surfaces.)


Figure 2: User controls and lights

**Relaxis** | [www.myrelaxis.com](http://www.myrelaxis.com)  
 1235 Puerta del Sol #600 | San Clemente, CA 92673 USA  
 T 1-888-475-7435 | F 1-949-492-0650

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